

Liquivida and

NIAGENTM Plus

Niagen IV is the pinnacle of pharmaceutical-grade NAD⁺ therapy, delivering you the best NAD⁺ boosting infusion experience for peak cellular health support.



Click here

Watch Webinar

*Cannot ship to California.

Niagen, the patented nicotinamide riboside chloride (NRCI), powers our next-generation treatment, supporting cellular energy, repair, and performance to support healthier aging.



Designed for those who prioritize efficiency, Niagen IV delivers superior results in a fraction of the time.



Each vial is \$195

It will be sold in a 4 pack.


Total \$780
Shipping included.

NiagenPlus 500 mg IV


ADMINISTRATION INSTRUCTIONS



Preparation



Reconstitute each vial of NiagenPlus with 5 mL of sterile water or from the Saline bag.



Transfer the reconstituted solution into a 500 cc saline bag.



Gently agitate the saline bag to ensure proper mixing and to prevent formation of a bolus.

Administration Options

Option 1 (Titrated Drip):

- Begin infusion at 50% rate.
- Increase speed gradually every 10 minutes until desired rate is achieved.

Option 2 (Rapid Infusion):

- Open infusion line fully ("wide open") and administer from start to finish.

Duration:

Total infusion time should be **20-40 minutes**, depending on the individual patient's tolerance and clinical response.

Note: It is solely the Practitioner's or Physician's responsibility to determine the appropriate administration, dosing, and frequency, as any guidelines provided are only suggestions.

Niagen Plus IV

Q&A GUIDE



What does each vial contain?

Each **10 mL vial** contains **500 mg** terminally sterilized by e-beam NRCL. **Refrigerated NAD+ (nicotinamide riboside chloride – NRCL).**

How often can it be administered?

Once per week.

Can it be paired with any other IV?

Yes – it can be combined with other IVs, **just not in the same bag as Glutathione.** Glutathione and Niagen are a **great combination**, but must be **run in separate bags** (drip one after the other).

What size saline bag should be used?

Either **500 mL or 1,000 mL** saline bag.

What is the recommended drip rate?

- **500 mL bag:** Minimum **3 drops/sec**, Maximum **5 drops/sec**
- **1,000 mL bag:** Minimum **6 drops/sec**, Maximum **10 drops/sec**

What is the retail value?

\$400-\$600 per IV treatment.

Can Niagen be made into a multi-dose vial (MDV)?

Yes. Reconstitute with **10 mL of Bacteriostatic Saline.**

Can Niagen be used for IM or Sub-Q injections?

Yes. Suggested dose is **0.5 mL (50 mg)** per injection. **Suggested retail: \$45 per injection.**

Can you use less dosing for IV administration?

Yes. For a lighter IV dose, use **5 mL from the MDV vial** for a total dose of **250 mg.**

Note: It is solely the Practitioner's or Physician's responsibility to determine the appropriate administration, dosing, and frequency, as any guidelines provided are only suggestions.